







Aplikasi :

- Failure Textile Inspection
- Composite Fiber Inspection
- Coloring
- Contaminate
- Thread Count Textile
- Etc ...

Mikroskop digital Dino-Lite memungkinkan para profesional industri Tekstil dan Komposit untuk melihat lebih banyak dan melakukan lebih banyak dalam menangkap berbagai gambar dan video dengan bantuan perangkat lunak.

Banyak profesional menggunakan Dino-Lite untuk melakukan inspeksi, tetapi ada berbagai model untuk dipilih dengan fitur yang berbeda tergantung pada kebutuhan pengguna.



High Series 200X



AM7915MZT
10~220X | 5MP | USB2.0
ESD Safe | Metal Housing
Measuring | Gambar & Video | EDOF | EDR



RK-10VXEXFPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)



MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Mid Series 950X



AM7515MT8A
700~900X | 5MP | USB2.0 | FLC
Measuring | Gambar & Video
Brightfield / Darkfield | AMR



RK-10VXEXFPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)



MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Basic Series 200X



AM4515ZT
20~220X | 1.3MP | USB2.0
ESD Safe | Plastic Housing
Measuring | Gambar & Video | Polarizer



RK-10VXEXFPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)



MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Eyepiece Series

untuk Trinocular



Euromex Calibration Scale
AE-110CALIB 1mm/100 (AE-1110)



Dino-Eye AM7025X
5MP | USB2.0
Measuring | Gambar & Video
23mm, 30mm, 30.5mm diameter adaptor