



Aplikasi :

- *Pes Control*
- *Wood Working*
- *Watch Repair*
- *Medical Transfusion*
- *Skin Inspection*
- *Etc ...*

Mikroskop Dino-Lite memberikan kemudahan para pengguna profesional pada fungsi pengukuran dengan gambar yang berkualitas tinggi.

Untuk item yang mungkin memiliki permukaan reflektif, kami merekomendasikan unit dengan pencahayaan terpolarisasi. Fitur polarisasi digunakan untuk mengurangi silau dan membantu menghilangkan 'titik panas' dari permukaan reflektif. Ini membantu memberikan pencahayaan seragam untuk membawa detail permukaan lebih banyak dengan silau minimal.

High Series 950X



AM7515MT8A
700~900X | SMP | USB2.0 | FLC
Measuring | Gambar & Video
Brightfield | Darkfield | AMR



RK-10VXEXFPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)



MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Mid Series 200X



AM7915MZT
10~220X | SMP | USB2.0
ESD Safe | Metal Housing
Measuring | Gambar & Video | ED OF | EDR



RK-10VXEXFPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)



MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Basic Series 200X



AM4515ZT
20~220X | 1.3MP | USB2.0
ESD Safe | Plastic Housing
Measuring | Gambar & Video | Polarizer



RK-10VXEXFPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)



MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Eyepiece Series

untuk Trinocular



Euromex Calibration Scale
AE-110CALIB 1mm/100 (AE-1110)



Dino-Eye AM7025X
5MP | USB2.0
Measuring | Gambar & Video
23mm, 30mm, 30.5mm diameter adaptor