






Kamera mikroskop Dino-Lite memungkinkan pengajar untuk menampilkan objek-objek yang kecil pada layar monitor yang besar. Ini dapat digunakan dalam pengaturan ruang kelas besar untuk kelompok besar siswa sebagai alat pengajaran kelompok yang bermanfaat.



Untuk tujuan pendidikan dan pengajaran, Dino-Lite adalah kamera mikroskop yang sempurna karena siswa dan guru dapat melakukan pemeriksaan mata pelajaran yang cermat dan berbagi dengan sekelompok besar rekan kerja atau teman sekelas.



- Aplikasi :**
- Breeding Spores Research
 - Zebra Fish Research
 - Nanoparticle Research
 - Stem Cell Research
 - Fluorescence Guide Surgery
 - Pest Control
 - Etc ...



High Series 200X

AM7915MZT
10~220X | 5MP | USB2.0
ESD Safe | Metal Housing | Measuring Gambar & Video | EDOF | EDR

RK-10VXEXFXPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)

MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Mid Series 950X

AM4515T8
700x~900x | 1.3MP | USB2.0
AMR | Plastic Housing | Measuring Gambar & Video | Magnification Locked

RK-10VXEXFXPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)

MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Basic Series 200X

AM4515ZT
20~220X | 1.3MP | USB2.0
ESD Safe | Plastic Housing
Measuring | Gambar & Video | Polarizer

RK-10VXEXFXPX
Professional Stand
22 cm (L) x 15 cm (W) x 23 cm (H)

MS15X-S1
Inspection Stage | 5 cm diameter
10 cm long x 8 cm wide | Vernier Scale

Eyepiece Series
untuk Trinocular

Euromex Calibration Scale
AE-110CALIB 1mm/100 (AE-1110)

Dino-Eye AM7025X
5MP | USB2.0
Measuring | Gambar & Video
23mm, 30mm, 30.5mm diameter adaptor